

# SOP: Competent cell preparation

By

**Richard Yu**



Research fellow



Centre for Marine Environmental Research and Innovative Technology (MERIT)  
Department of Biology and Chemistry  
City University of Hong Kong  
(Email: [bhrichyu@cityu.edu.hk](mailto:bhrichyu@cityu.edu.hk))

## **Day 1**

### ***Procedure:***

1. Streak out the DH5 $\alpha$  frozen stock (stored in  $-80^{\circ}\text{C}$  freezer) onto a LB plate (without antibiotics). Incubate overnight (16–20 hours) at  $37^{\circ}\text{C}$ .
2. Autoclave 200 ml (in 250-ml blue-cap bottle) and 500 ml (in a 1-L Erlenmeyer flask)  $\psi$  broth for Day 2 experiment.

## **Day 2**

### ***Procedure:***

1. Transfer 25 ml autoclaved  $\psi$  broth (from the 250-ml blue-cap bottle) to a sterile 100-ml Erlenmeyer flask and inoculate with a single DH5 $\alpha$  colony on the plate. Grow overnight (no more than 16 hours) in a  $37^{\circ}\text{C}$  shaker (250 rpm).
2. For Day 3 experiment, place the 1-L flask containing 500 ml autoclaved  $\psi$  broth in a  $37^{\circ}\text{C}$  incubator.

## **Day 3**

### ***Things to prepare:***

1. Switch on the UV-VIS spectrophotometer (Shimadzu UV-1206) and warm it up for at least 30 min.
2. Switch on the Beckman Avanti J-25I centrifuge and maintain it at  $4^{\circ}\text{C}$ .
3. Chill two 500-ml centrifuge bottles on ice.
4. Pre-cool 200 pcs autoclaved 1.5-ml eppendorf tubes.
5. Get 1–2 liters of liquid nitrogen.

***Procedure:***

1. Take out 1 ml  $\psi$  broth from the 1-L flask and save it (the blank for OD measurement).
2. Transfer the 25 ml overnight culture to the 1-L flask. Swirl the flask to mix and measure the initial OD<sub>550</sub> of the culture.
3. Incubate the culture in a 37°C shaker (250 rpm).
4. Assume that the bacteria will double every 20 min, estimate the shortest time required for the culture to reach OD<sub>550</sub> = 0.44 (it takes approx. 45 min to 1 h). Thereafter, measure OD<sub>550</sub> every 10–20 min.
5. When the culture OD<sub>550</sub> reaches 0.44, transfer the flask on an ice-bath (Note: It is important to keep the cells at 4°C (or on ice) for the remainder of the procedure).
6. Split the culture into two equal parts by pouring ~250 ml of the culture into each chilled 500-ml centrifuge bottle. Measure and balance the weights of the bottles using a digital balance.
7. Pellet the cells at 5,000 rpm at 4°C for 10 min using the Beckman Avanti J-25I centrifuge.
8. Decant the supernatant (into the emptied 1-L Erlenmeyer flask) and resuspend thoroughly the cell pellet (on ice) in 2 ml chilled TfbI using a 1000- $\mu$ l filter tip. Add additional 98 ml TfbI to each bottle and mix by gentle swirling.
9. Hold the content on ice for 5 min.
10. Pour the contents together. Centrifuge the cells at 4,000 rpm at 4°C for 10 min.
11. Decant the supernatant and resuspend thoroughly the cells in 4 ml chilled TfbII using a 1000- $\mu$ l filter tip. Add additional 16 ml TfbII to the cells and mix by gentle swirling.
12. Hold the content on ice for 15 min.
13. Aliquot 100  $\mu$ l (use a 1000- $\mu$ l tip) into each pre-cooled 1.5-ml eppendorf tube.
14. Quickly freeze the cells in nitrogen liquid prior to storage in a –80°C.

**Determination of Transformation Efficiency of Competent Cells**

***Things to prepare:***

1. Set a water bath at 42°C.
2. Take out 4 LB agar/amp plates from 4°C refrigerator and prewarm in a 37°C air-incubator.

***Procedure:***

1. Serially dilute a plasmid [e.g. pGEM 3Z(f+)] stock to 20 pg/ $\mu$ l and 2 pg/ $\mu$ l. Below is an example for such a serial dilution:

$$\text{pGEM 3Z(f+)} \text{ stock} = 0.2 \text{ g/L} = 0.2 \text{ } \mu\text{g}/\mu\text{l} = 200 \text{ ng}/\mu\text{l} = 200 \text{ 000 pg}/\mu\text{l}$$

***Dilution A:***

1  $\mu$ l pGEM 3Z(f+) stock (use a P2.5 autopipette and a P10 tip to dispense)  
+ 99  $\mu$ l water  $\rightarrow$  2000 pg/ $\mu$ l

***Dilution B:***

10  $\mu$ l Dilution A + 90  $\mu$ l water  $\rightarrow$  200 pg/ $\mu$ l

***Dilution C:***

10  $\mu$ l Dilution B + 90  $\mu$ l water  $\rightarrow$  **20 pg/ $\mu$ l**

***Dilution D:***

10  $\mu$ l Dilution C + 90  $\mu$ l water  $\rightarrow$  **2 pg/ $\mu$ l**

2. Remove 4 tubes of frozen competent cells from  $-70^{\circ}\text{C}$  freezer and place them on ice until just thawed (about 5 mins). Label them as **100 pg (Replicate 1)**, **100 pg (Replicate 2)**, **10 pg (Replicate 1)** and **10 pg (Replicate 2)**.
3. Add 5  $\mu$ l **Dilution C** to the 100 pg (Replicate 1) and 100 pg (Replicate 2) tubes and 5  $\mu$ l **Dilution D** to the 10 pg (Replicate 1) and 10 pg (Replicate 2).
4. Gently flick the tube, or use a pipette tip to stir, to mix and place the tubes on ice for 20 minutes. **AVOID PIPETTING UP AND DOWN TO MIX, AS THE CELLS ARE VERY VERY FRAGILE.**
5. Heat shock the cells for 50 seconds in a water bath at exactly  $42^{\circ}\text{C}$ .
6. Immediately return the tubes to ice for 2 minutes.
7. Add 1000  $\mu$ l of room temperature SOC or LB medium to the tubes.
8. Incubate for 1.5 hours in a bath water at  $37^{\circ}\text{C}$  with or without shaking (ca. 150 rpm).
9. Centrifuge the cells at 3000 rpm for 10 mins using a microcentrifuge. Decant the supernatant and resuspend the cells in about 50  $\mu$ l LB broth.
10. Plate them onto LB agar/amp plates using a hockey stick.
11. Incubate the plates overnight (16 hours) at  $37^{\circ}\text{C}$ . Calculate the average transformation efficiency based on the colony counts on the plates. This should be at  $1 \times 10^6 - 1 \times 10^7$  cfu/ $\mu$ g plasmid DNA (i.e. 100–1000 cfu on the 100 pg plates and 10–100 cfu on the 10 pg plates)

## Appendix

### $\psi$ broth:

	<u>500 ml</u>	<u>200 ml</u>	<u>Final conc.</u>
Tryptone	10g	4 g	2%
Yeast extract	2.6 g	1.04 g	0.5%
MgSO <sub>4</sub>	2.4 g	0.96 g	20 mM
NaCl	0.3 g	0.12 g	10 mM
KCl	0.94	0.376 g	25 mM

Autoclave at 121°C for 15 min.

### TfbI buffer (per 200 ml):

	<u>Amount</u>	<u>Final conc.</u>
<b>KAc</b>	<b>0.59 g</b>	<b>30 mM</b>
<b>CaCl<sub>2</sub></b>	<b>0.29 g</b>	<b>10 mM</b>
(1) Stir to dissolve and adjust pH to 5.8 with acetic acid (only a few drops).		
<b>MnCl<sub>2</sub></b>	<b>1.98 g</b>	<b>50 mM</b>
<b>RbCl<sub>2</sub></b>	<b>2.42 g</b>	<b>100 mM</b>

(2) Stir to dissolve.

**Glycerol 30 ml 15%**

(3) Make the final volume to 200 ml with deionized H<sub>2</sub>O.

(4) Filter sterile the buffer using a 50-ml syringe attached to a 0.45- $\mu$ m syringe filter.

(5) Collect the filtrate with a sterile 250-ml CC bottle

(6) Store the buffer at 4°C.

### TfbII buffer (per 200 ml):

	<u>Amount</u>	<u>Final conc.</u>
<b>MOPS</b>	<b>0.462 g</b>	<b>10 mM</b>
<b>CaCl<sub>2</sub></b>	<b>2.36 g</b>	<b>75 mM</b>
<b>RbCl<sub>2</sub></b>	<b>0.24 g</b>	<b>10 mM</b>

(1) Stir to dissolve and adjust pH to 6.5 with KOH (only a few drops).

**Glycerol 30 ml 15%**

(3) Make the final volume to 200 ml with deionized H<sub>2</sub>O.

(4) Filter sterile the buffer using a 50-ml syringe attached to a 0.45- $\mu$ m syringe filter.

(5) Collect the filtrate with a sterile 250-ml CC bottle

(6) Store the buffer at 4°C.