

SOP: RNA Blot Hybridization

By

Richard Yu



Research fellow



Centre for Marine Environmental Research and Innovative Technology (MERIT)
Department of Biology and Chemistry
City University of Hong Kong
(Email: bhrichyu@cityu.edu.hk)

A. Electrophoresis of RNA

1. Preparative work:

- i) Clean the followings with laboratory detergent and soak in 3% hydrogen peroxide (H_2O_2) for 20 min:
 - electrophoresis apparatus
 - Erlenmeyer flask
 - measuring cylinders (50, 500, 1000 ml)
 - plastic container
 - blot transfer tray

Rinse thoroughly with Milli-Q H_2O (or DEPC H_2O) and dry.

- ii) Deionize formaldehyde and formamide. It is a must when formaldehyde looks yellow in color or has a $\text{pH} < 4$ (i.e. oxidized). Formamide in the deionized form can be available from Amresco. *Deionization*: Add 1 g of mixed bed resin (AG 501-X8, Bio-Rad or M-8032, Sigma) to 10 ml of formaldehyde/formamide. Swirl the mixture every 3-5 min within the 1 h-deionization. Add extra resin when the beads turn golden-yellow (i.e. an indication of saturation). Remove the resin by filtration. The deionized materials should be used immediately. Discard any unused portion.

- iii) Prewarm appropriate amounts (see below) of the 10×MOPS buffer (0.2 M MOPS, pH 7; 50 mM NaOAc; 10 mM EDTA, pH 8.0; filter sterilize and store in a brown bottle) and 36.5 % formaldehyde in a 55-60°C waterbath. Reserve a heat-block or PCR machine if possible. Empty a place in the fume hood for pouring the gel.
2. Prepare 200 ml (50 ml) of molten agarose solution by melting 2.4 g (0.6 g) of agarose (use a non-foaming one, such as BioRad or Promega) in 170 ml (42.5 ml) DEPC H₂O in an Erlenmeyer flask in a microwave oven. Note: After adding agarose to water, wait for at least 1 min before melting (this can minimize foaming). Allow to cool to 55-60°C under running water and place in the preset waterbath (at 55-60°C). Equilibrate for 15 mins.
 3. Add 20 ml (5 ml) prewarmed 10×MOPS buffer and 10 ml (2.5 ml) prewarmed 36.5% formaldehyde. This will produce a solution of 1.2% agarose in 1×MOPS buffer and 0.6 M formaldehyde.
 4. Cast a 0.5- to 0.75-cm-thick gel (a thinner gel offers better resolution of the RNA) using the solution prepared in (3). This should be done in a fume-hood to minimize formaldehyde fumes exposure (formaldehyde is toxic and teratogenic). Allow the gel to solidify for at least 30 min.
 5. Prepare a 1.25×RNA loading buffer as follows:

100% formamide	62.5 μl
10×MOPS	12.5 μl
36.5% formaldehyde	24.0 μl
<u>4% bromophenol blue</u>	<u>1.0 μl</u>
Total:	100 μl (for 10 RNA samples)

6. Prepare the RNA for denaturation by mixing the following in a sterile 0.2 or 0.5 ml thin-wall microcentrifuge tube:

2.5 μl RNA (up to 20 μg; dissolved in DEPC H₂O)
 10.0 μl RNA loading buffer
 0.5 μl ethidium bromide (1 mg/ml; prepared with DEPC H₂O)

Thorough pipetting is essential.

Note:

- i) Do not exceed 20 μg RNA per lane to prevent overloading and losing resolution. If rare mRNAs are of interest, as much as 4-5 μg of poly(A)⁺ RNA can be applied per lane. But in any case, use of mRNA (e.g. 1 μg for moderate expression) always produces better resolution and lower background.

- ii) An appropriate RNA marker (e.g. Gibco/BRL 0.16-1.77 kb or 0.25-0.95 kb RNA ladder) should be included.
7. Denature all RNA samples and the RNA marker by heating at 65°C for 10 min using a PCR thermal cycler. Quench on ice for 2 min. Quick spin to collect the mixture at the bottom of the tube.

During the incubation of the RNA, immerse the gel in running buffer (1×MOPS; without ethidium bromide). Note: Use the minimal amount of buffer necessary to completely cover the gel.

8. Electrophorese the samples at a maximum of 5 V/cm distance between the electrodes (usually 100 V is OK). At the end of electrophoresis, capture the image with inclusion of dimension (using a fluorescent ruler) under an UV transilluminator.
9. Soak the gel in plenty of DEPC-H₂O (or Milli-Q H₂O) and apply gentle shaking to remove formaldehyde. Typically, a 30 min-soaking in several changes of H₂O is adequate.

B. Northern Transfer

10. Capillary transfer can be a choice for RNA blotting. Transfer (5×SSC) should be allowed to take place overnight (i.e. at least 16 h). Baking at 80°C for 2 h is a good means for immobilization and storage (either at room temp. and away from moisture or at -20°C). Once the RNA has been immobilized, it is relatively resistant to degradation.

C. Hybridization

11. Hybridization using ExpressHyb (Clontech) is highly recommended. Hybridization can be performed at 60-68 °C for 2 h using a cDNA probe at 2×10⁶ cpm/ml. Membrane is washed according to the manufacturer's instructions (i.e. initial wash: 2×SSC, 0.05% SDS, rm temp for 30-40 min with several changes; final wash: 0.1×SSC, 0.1% SDS, 50°C for 40 min with one change).

D. Autoradiography

12. Expose the membrane to the K-screen (Kodak) for overnight and capture the image using the Fx molecular imager (BioRad). Over-exposure should be prevented if quantification is required (the Fx software will indicate saturated pixels).
13. For normalization, strip the membrane accordingly and probe with a housekeeping gene (e.g. β-actin) or a rRNA gene.