

SOP: Primer Extension

By

Richard Yu



Research fellow



Centre for Marine Environmental Research and Innovative Technology (MERIT)
Department of Biology and Chemistry
City University of Hong Kong
(Email: bhrichyu@cityu.edu.hk)

A. Oligonucleotide labeling by phosphorylation

1. Add the following components to a 1.5-ml microcentrifuge tube:

Oligonucleotide primer (2 μ M)	5 μ l (10 pmol)
5 \times Forward Reaction Buffer	10 μ l
T4 Polynucleotide Kinase (Gibco BRL)	2 μ l (20 units)
[γ - ³² P]ATP (10 μ Ci/ μ l, 3000 Ci/mmol)	5 μ l (50 μ Ci)
Sterile water	28 μ l
	Total: 50 μ l

Note:

- a) The oligonucleotide used should be:
 - (i) a 30- to 40-mer to allow a reasonable stringency of hybridization conditions; and
 - (ii) exon I sequence, close to the predicted transcription start site (not more than 100–150 bp).
 - b) Instead of [α -³²P]dCTP, [γ -³²P]ATP must be used in this analysis.
2. Incubate at 37°C for 10 min.
 3. Heat at 65°C for 10 min to stop the reaction.
 4. Purify the oligonucleotide using the QIAquick Nucleotide Removal Kit from Qiagen. Elute the oligonucleotide with 200 μ l EB buffer.

5. Count 5 μl by liquid scintillation counting.

B. Hybridization

1. Mix 2×10^4 cpm of the labeled oligonucleotide with 10 μg (up to 50 μg , depending on gene expression level) of the total RNA analyzed (Note: Overloading of primer would cause artifactual bands).
2. Add 0.1 volume of 3 M sodium acetate (pH 5.2) and 2.5 volumes of absolute ethanol (-20°C). Store the solution at -80°C for at least 1 h.
3. Centrifuge at $12,000 \times g$ for 15 min at 4°C .
4. Wash the pellet with 70% ethanol (-20°C) and re-centrifuge for 5 min at 4°C .
5. Carefully remove all of ethanol and air-dry the pellet for 5 min (Note: An over-dried pellet cannot be easily re-dissolved).
6. Redissolve the pellet in 30 μl of hybridization buffer (40 mM PIPES (pH 6.4); 1 mM EDTA (pH 8.0); 0.4 M NaCl; 80% deionized formamide) by drawing through a pipettor tip 50–100 times.
7. Denature the nucleic acids at 85°C for 10 min. Quickly transfer the hybridization mixture to a water bath set at $30\text{--}35^\circ\text{C}$. Incubate for 8–12 h or overnight (i.e. 16 h).
8. Add 170 μl sterile water and 400 μl absolute ethanol (-20°C), mix well, and place on ice for 1 h. Repeat steps 3–5.

C. Primer extension

1. Redissolve the RNA:primer hybrid in the following RT buffer:

5 \times First Strand Buffer	8 μl
0.1 M DTT	4 μl
50 mM dNTP mix	4 μl
RNase inhibitor or RNaseOut (Gibco, BRL)	2 μl
DEPC-treated water	20 μl

Note: Such a high dNTP concentration (i.e. 5 mM) helps minimize artifact bands.

2. Incubate at 42°C for 2 min. Add 2 μl (400 U) of SUPERSCRIT II and gently re-suspend. Incubate for 50 min at 42°C .

3. Inactivate the reaction by heating at 70°C for 15 min. Slowly cool down to room temperature.
4. Add 2 µl DNAase-free pancreatic RNAase (5 µg/ml). Incubate the reaction for 30 min at 37°C. A further 3 min-incubation at 95°C can enhance complete RNA degradation. Slowly cool down to room temperature.
6. Add 158 µl TE buffer and 200 µl phenol: chloroform:IAA (24:24:1). Vortex for 30 s, and centrifuge at full speed for 5 min at room temperature.
6. Transfer the upper, aqueous phase to a clean tube. Add 0.1 volume of 3 M sodium acetate (pH 5.2), 2.5 volumes of absolute ethanol (−20°C) and glycogen.
7. Store the solution at −80°C for at least 1 h. Centrifuge at 14,000 rpm for 15 min at 4°C. Wash the pellet with 70% ethanol (−20° C) and re-centrifuge for 5 min at 4°C. Air-dry the primer-extended product for 5 min.
8. Dissolve the pellet in 4 µl TE buffer, add 6 µl formamide loading buffer (80% formamide, 10 mM EDTA (pH 8.0), 1 mg/ml xylene cyanol FF, 1 mg/ml bromophenol blue), and mix well. Store at −20°C before gel electrophoresis.
9. Heat 5 µl of the mixture for 5 min at 95°C. Immediately chill on ice.
10. Analyze the radiolabeled cDNA by co-electrophoresis with a DNA sequencing ladder generated from an appropriate genomic clone (exon I + promoter) using the same oligonucleotide of above, through a 6% denaturing polyacrylamide gel.