

SOP: SYBR Green-based real-time RT-PCR

By

Richard Yu



Research fellow



Centre for Marine Environmental Research and Innovative Technology (MERIT)
Department of Biology and Chemistry
City University of Hong Kong
(Email: bhrichyu@cityu.edu.hk)

A. Development and assessment of assays

1. **Design of gene-specific primers**

- **Primer3**

http://frodo.wi.mit.edu/cgi-bin/primer3/primer3_www.cgi

This is used for picking up primers from an input gene sequence. Below are the recommended settings for selecting real-time PCR primers:

- “*Product size range*” = 80–150 bp
- “*Primer size*” = Min: 18 Opt: 20 Max: 23
- “*Primer Tm*” = Min: 58 Opt: 60 Max: 63
- From the output primer list, select the primers which end with G/C at their 3'-ends or avoid the primers ended with several consecutive A/T at their 3'-ends.

- **Oligoanalyzer**

<http://www.idtdna.com/analyzer/Applications/OligoAnalyzer/Default.aspx>

This is used for more detailed analyses for your selected primers for undesirable structures and specificity (e.g. hairpin, self-dimer, hetero-dimer and BLAST check).

2. **Validation of primer efficiency and specificity by dilution curve construction**

- Prepare 4 concentrations of tissue cDNA template by serially diluting a standard RT reaction product (see Appendix 1) or plasmid DNA (10^4 – 10^7 copies/reaction) as follows:

<u>Dilution</u>	<u>Recipe</u>
A (1/5):	4.5 μ l RT rxn product + 18 μ l water
B (1/20):	6 μ l Dilution A + 18 μ l water
C (1/80):	6 μ l Dilution B + 18 μ l water
D (1/320):	6 μ l Dilution C + 18 μ l water

Add 5 μ l of each concentration in triplicate to appropriate wells on a 96-well PCR plate.

	Triplicate		
Conc.	1	2	3
A	○	○	○
B	○	○	○
C	○	○	○
D	○	○	○

- Prepare a PCR master mix as follows:

	<u>1 rxn. (μl)</u>	<u>Master mix 13 rxn. (μl)</u>
2 \times iQ SYBR Green Supermix*	12.5	162.5
10 μ M Forward primer	0.5	6.5
10 μ M Reverse primer	0.5	6.5
Water	6.5	84.5
Diluted cDNA (A, B, C, or D)	5	(pre-add to well)
	Total = 25	

* For other PCR machines which require ROX for normalization (e.g. ABI 7500), use other ROX-containing supermixes (e.g. ABI SYBR Green Master Mix).

3. Mix 20 μ l PCR master mix with each cDNA aliquot by gently pipetting up and down several times (Note: Try your best to minimize bubble formation).
4. Seal the 96-well plate with an optical tape (BioRad, Cat# 2239444 or equivalent). Spin down all droplets by centrifugation at 4,000 rpm for 2 min.


5. Warm up the iCycler system for 15–20 min before use. Use the following protocol for a general SYBR Green PCR assay:

Cycle 1: (1X)
 Step 1: 95.0°C for 00:30
Cycle 2: (45X)
 Step 1: 95.0°C for 00:20
 Step 2: 60.0°C for 00:30
 Step 3: 72.0°C for 00:30
 Data collection and real-time analysis enabled.

Melt curve analysis:
Cycle 3: (1X)
 Step 1: 95.0°C for 01:00
Cycle 4: (1X)
 Step 1: 55.0°C for 01:00
Cycle 5: (80X)
 Step 1: 55.0°C for 00:10
 Increase setpoint temperature after cycle 2 by 0.5°C
 Melt curve data collection and analysis enabled.

6. **Typical results of a successful assay:**

a. PCR amplification plot (PCR Base Line Subtracted)




Note: The threshold line should be placed within the **exponential increase phase** (which looks linear in the log transformation) of the amplification curve. Check from the plot that the curves **DO NOT** overlap with the **baseline range** (e.g. preset to be **2–10 cycles** in iCycler).

b. Standard curve (PCR Base Line Subtracted)

Correlation Coefficient: 0.999 Slope: -3.388 Intercept: 22.287 $Y = -3.388 X + 22.287$
PCR Efficiency: 97.3 %


□ Unknowns
○ Standards



PCR efficiency and R^2 are calculated automatically by the iCycler software from the dilution curve. Otherwise, use the following formula for the calculation of PCR efficiency in other real-time PCR systems:

$$\text{PCR efficiency} = (10^{-1/\text{slope}} - 1) \times 100\%$$

c. Melt-curve analysis



Specific amplification should be represented by a **single peak** emerged at about 80–90°C. Minor peaks at lower and higher temperatures implies the presence of primer dimers and larger non-specific PCR products, respectively.

7. Optimization of PCR reaction

- A. *If PCR efficiency is less than 90% (i.e. low sensitivity of the assay):*
 - Increase primer concentration from 200 nM to 400–600 nM; or
 - Re-design primers
- B. *If more than 1 peak in the melt curve:*
 - Decrease primer concentration from 200 nM to 100 nM if primer dimer formation occurs; or
 - Increase annealing temperature from 60°C to 62°C if other non-specific product form; or
 - Re-design primers

C. *If target gene expression in the selected tissue is too low for dilution curve construction:*

- Use plasmid DNA containing the target gene insert as a template (concentration range: 10^6 , 10^5 , 10^4 and 10^3 copies/rxn) for dilution curve construction.
- Use the following formula for copy number calculation:
No. of molecules/ μ l
 $= (\text{plasmid conc. in ng} \times 1.515 \div N_{\text{bp}}) \times 6.023 \times 10^{11}$
where N_{bp} = size of vector + size of insert in bp

8. **Verification of amplification specificity**

- Agarose gel electrophoresis (there should be one band seen on the gel)
- Cloning of the PCR product for automated DNA sequencing
Note: Don't clone PCR product generated from UDG-containing PCR supermixes because bacteria fail to recognize using uracil-containing template for DNA replication. If you use UDG-containing PCR supermixes, direct sequencing of your PCR product is the only choice.

B. Gene expression measurements


1. Dilute RT product (i.e. cDNA) by 10 times as follows:
1.6 μ l RT product + 14.4 μ l H₂O

If the expression of your target gene is low, dilute RT-product by 5 times. In general, 1/5 dilution is the highest concentration which you can input for PCR. Concentrations higher than 1/5 will not be linearly related to C_T as predicted from the dilution curve.

2. Add 5 μ l of the diluted RT product in triplicate to the wells of a 96-well PCR plate.
3. See above for other procedures and record the C_T values for samples (Use the same threshold for all samples).

C. 18S rRNA Normalization

1. Dilute RT products by 5000 times as follows:
 - (a) Prepare a 1/100 dilution by adding 1 μ l RT product (dispense using a P20 autopipette) to 99 μ l water. Mix thoroughly by vortexing.
 - (b) Transfer 2 μ l of the diluted RT product to 98 μ l water aliquoted in a new tube. Mix thoroughly by vortexing. This tube contains RT product at 1/5000 dilution.



2. Add 5 µl of the diluted products to the wells of a 96-well PCR plate.
3. Prepare a master mix as follows:

	<u>1 rxn. (µl)</u>	Master mix <u>x+1 rxn. (µl)</u>
iQ SYBR Green Supermix	12.5	
10 µM om18S-F	1	
10 µM om18S-R	1	
Water	5.5	
1/5000 diluted cDNA	5	(pre-add to well)
	Total = 25	

Oryzias melastigma-specific 18S rRNA primer:
om18S-F (5'-CCTGCGGCTTAATTTGACCC-3')
om18S-R (5'-GACAAATCGCTCCACCAACT-3')

4. Add 20 µl of master mix to each well. Mix by pipetting up and down for several times.
5. See above for other procedures and record the C_T values for samples (Use the same threshold for all samples).

D. Calculation of gene expression

For each experimental sample, the expression level of a target gene was normalized to 18S rRNA, resulting in a mean normalized expression (MNE) value as derived from Equation (1):

$$MNE = \frac{(E_{\text{reference}})^{C_{T\text{reference,mean}}}}{(E_{\text{target}})^{C_{T\text{target,mean}}}} \quad (1)$$

where E_{target} and $E_{\text{reference}}$ ($=10^{-1/\text{slope}}$) (Simon, 2003). Final data were calculated as mean normalized target gene expression relative to the control (Equation 2):

$$\text{N-fold change} = \text{MNE}_{\text{exp}}/\text{MNE}_{\text{con}} \quad (2)$$

Reference

Simon, P., 2003. Q-Gene: processing quantitative real-time RT-PCR data. *Bioinformatics* 19, 1439–1440.

Appendix 1 First-strand cDNA synthesis by reverse transcription

A. DNase Digestion

1. Set up the DNase digestion as follows:

Total RNA	1 μg (in 8 μl H ₂ O)
10x Reaction buffer*	1 μl
DNase (1 U/ μl)*	1 μl
	Total: 10 μl

*provided in the RQ1 RNase-Free DNase kit (Promega, Cat.# M6101)

2. Incubate at 37°C for 30 min.
3. Add 1 μl Stop Solution to terminate the reaction.
4. Incubate at 65°C for 10 min to inactivate the DNase.

B. Reverse transcription

1. Add 3 μl random primer (40 ng/ μl) to each tube containing the DNase digestion product.
2. Heat to 70°C for 5 min. Cool quickly on ice for 5 min.
3. Prepare a “**RT master mix**” as follows:

	1 reaction	Master mix <i>n</i> +1 reaction
5x Reaction buffer*	5	
10 mM dNTP	1.25	
RNaseOut (Invitrogen)	1	
Water	2.75	
M-MLV RT (H-)*	1	

*M-MLV Reverse Transcriptase, RNase H Minus, Point Mutant (Promega, Cat.# M3683)

4. Add 11 μ l RT master mix to each tube in (3). Pipette up and down to mix the contents (total vol. = 25 μ l).
 5. Incubate at room temperature for 10 min.
 6. Incubate at 42°C for 50 min.
 7. Inactivate the reaction at 70°C for 15 min.
 8. Store at -20°C until use.
-

For more information, please go to my website, Molbio-Spot (<http://www.geocities.com/richyuhk/>) and see under PCR technology: Real-time PCR.