


# Top-down Analysis of Modifications on Cysteine 34 of Albumin in Human Serum Samples

Jianzhong Chen<sup>1</sup>, Nathan K Kaiser<sup>1</sup>, Hollie Huff<sup>2</sup>,  
Yvonne Carella<sup>2</sup>, James E Bruce<sup>1</sup>

<sup>1</sup>Washington State University, Pullman, WA

<sup>2</sup>Inverness Medical Innovations, Inc., Louisville, CO

# Human Serum Albumin (HSA)



- most abundant protein in serum
- 585 amino acid residues
- 35 cysteine residues:
  - 17 disulfide bonds
  - + 1 free cys34 with –SH group


Generated with Protein Explorer

[http://www.umass.edu/microbio/chime/pe\\_beta/pe/protexpl/frntdoo2.htm](http://www.umass.edu/microbio/chime/pe_beta/pe/protexpl/frntdoo2.htm)

## Modifications of –SH on Cys34 of HSA

- Thiols and disulfides in serum
  - cysteine, homocysteine, cysteinylglycine and glutathione
  - disulfides formed by thiols
- Cys34 modified by thiols or disulfides
  - maybe related to ischemia (restriction in blood supply )
  - potential biomarker for acute coronary syndromes
- Study Method
  - bottom-up & multiple reaction monitoring (MRM)

# Artificial Modification of Cys34 with Bottom-up Approach




Mock serum digestion experiment:


600  $\mu\text{M}$  synthetic  
Cys34 tryptic peptide (Cys34P)

6 mM Cysteine    50 mM  $\text{NH}_4\text{HCO}_3$   
37°C

Cys-Cys34P + Cys34P-Cys34P (Dimer)  
(Predominant)



# History of Top-down Analysis of Serum Albumin<sup>1,2</sup>



1. Smith RD et al, Anal Chem, 1990, 62, 882
2. Loo JA et al, Anal Chem, 1991, 63, 2488

# Top-down Approach

## Advantages:


- Direct Analysis
  - fast
  - decreased artificial modifications
- 100% sequence coverage
- Easier to identify protein isoforms and mutations

## Challenges:



- Low sensitivity
- Difficult fragment assignments
  - overlapping peaks
  - internal fragments
  - frequency shift due to space-charge effect

# Experimental


1  $\mu\text{L}$  Serum or 500  $\mu\text{M}$  HSA + 500  $\mu\text{L}$  spray solution  
(50:50:2 MeOH: H<sub>2</sub>O: HAc)



# MS Spectrum of HSA in Diluted Serum Sample




# Isotopic Resolution of Cys-HSA<sup>40+</sup> in Sigma HSA





High resolution detection makes it possible to do top-down analysis.

# MS/MS Spectrum of Sigma HSA<sup>44+</sup>



# Two Complementary Fragments of Sigma HSA<sup>44+</sup>



Monoisotopic  $M_r$  of Cys-HSA


Experimental: 66513.40 Da

Theoretical: 66513.66 Da


$\Delta M_r = -0.26$  Da

- 4.4 ppm


# Abundances of Different Forms of HSA in Serum Samples



# Cys34 Related HSA Fragments in Serum Samples



# Intensity Ratios of HSA Fragments in Serum Samples



## Summary

- Top-down approach is a more reliable method for the quantitation of modifications on Cys34 of HSA.
- Isotopic resolution and high mass accuracy MS/MS analysis of HSA were achieved for large fragments.
- Relative intensities of the fragments for different samples help assign the fragments or confirm the assignments.
- Modified and/or unmodified fragments may be used for relative quantitation.

# Acknowledgements

## Bruce Group at WSU

Professor James Bruce

Dr. Xiaoting Tang

Dr. Haizhen Zhang

Dr. Nathan Kaiser

Natalia Zakharova

Gunnar Skulason

Li Yang

Hye In Nam

Jake Zheng

Brian Webb

Chad Weisbrod

Megan Wolff

## Laboratory for Biotechnology and Bioanalysis at Washington State University

Dr. William Siems

Dr. Gerhard Munske

## Pacific Northwest National Lab

Gordon Anderson

## Inverness Medical Innovations, Inc

Dr. Hollie Huff

Dr. Yvonne Carella

## Funding

Inverness Medical Innovations, Inc.