



Research Scientific Services

MS-QuietBoxII noise reduction system

Tired of all that noise?
*Sound reduction
enclosure for your
vacuum pumps*



New MS-QuietBoxII

Improved solution to mass spec pump noise reduction!

- Reduce pump noise by 15 dB(A)
- Install without mass spec downtime
- 3 options for vacuum line access points for easy setup
- Integral fans maintains vacuum pump at optimum temperature (safety options available)
- Compatible with all common LCMS vacuum pumps
- Easy access for pump maintenance
- Flat pack delivery; can be assembled in 5 minutes - no tools required!
- Single or dual pump box available.

MS-QuietBoxI I

Ordering Information

| Pump make/type | Description | Part No. |
|---|--|-----------|
| Edwards RV3, RV5, XDS5, E2M0.7, Ulvac G100DC, Ulvac GLD040, Varian MS40+, Leybold S1, 5, D2, 5E and Triscoll 300/600 | MS-QuietBox II for pumps shown including tray | MSQB2-01T |
| Edwards RV8, RV12, E1M18, E2M18, E2M28, E2M30, Leybold Trivac D25B, D16B, D4B, D8B, D16B, Varian DS42, DS102, DS202, DS302, DS402, DS602, HS602, HS452, HS652, Alcatel 1005SD, 1015SD, 2004A, 2005SD/I/C1, 2010SD/I/C1, 2015SD/I/C1, 2021SD/I/C1, 2005H1, 2015H1, 2010C2, 2015C2, 202 1C2, 2002I & Duo 10 Pfeifer | MS-QuietBox II for pumps shown including tray | MSQB2-02T |
| Edwards XDS35I | MS-QuietBox II for Edwards XDS35I scroll pump | MSQB2-03 |
| Sogevac SV40B, SV60 | MS-QuietBox II for Sogevac SV40 including tray | MSQB2-03T |
| MS-QuietBox II Double Box | MS-QuietBox II double box for 2 pumps including tray | MSQB2-02D |
| Edwards E2M1.5, Duo 2.5 Pfeiffer | MS-QuietBox II for pumps shown | MSQB2-07 |



Research Scientific Services
PO Box 9188, Gaithersburg, MD 20989
301-977-9344

www.resci.com

RSS recognizes and acknowledges all trademarks. RSS has no OEM affiliation, E&OE.